

Fig. Schematic diagram of a spectrometer




Fig. Ray diagram of angle of prism

TOPIC:

Integrated Optics [Angle of Prism]

AIM:

To determine the apex angle of the given prism using a spectrometer.

APPARATUS:


- Spectrometer
- Spirit level
- Magnifying glass
- Glass prism
- Sodium vapour lamp

TABULATION:

$$\text{Least count} = \frac{0.5^\circ}{30 \text{ no. of div.}} = \frac{30'}{30} = 1'$$

| Reading of reflected ray | Vernier A | | | Vernier B | | |
|-----------------------------|-----------|-----|---------|-----------|-----|--------|
| | MSR | VSR | Total | MSR | VSR | Total |
| Reflection from side (a) AB | 114° | 14' | 114°14' | 324° | 8' | 324°8' |
| Reflection from side (b) AC | 263° | 9' | 263°9' | 263° | 9' | 263°9' |
| Difference between a & b | 118°55' | | | 118°1' | | |

$$\text{Mean } \theta = 59^\circ 14'$$



Umapati
08/08/19
09BCE2105

RESULT:

The Apex angle of the given equilateral prism = 59°14'

Teacher's Signature _____