Construction of Daniel cell and measurement of its voltage with three different concentrations of Cu/Zn solutions:

In the Daniel cell, copper and zinc electrodes are immersed in the equimolar solution of CuSO<sub>4</sub> and ZnSO<sub>4</sub> respectively.

At the anode, zinc is oxidized as per the following halfreaction:  $Zn_{(s)} \rightarrow Zn^{2+}_{(aq)} + 2e^{-}$ 

At the cathode, copper is reduced as per the following reaction:  $Cu^{2+}_{(aq)} + 2e^{-} \rightarrow Cu_{(s)}$ 

The overall reaction is:  $Zn_{(s)} + Cu^{2+}_{(sq)} \rightarrow Zn^{2+}_{(sq)} + Cu_{(s)}$ Construct Daniel cell using the following concentrations of Copper and Zinc solutions and record the voltage of the cells in Table 3.




Table 3: EMF of Daniel Cell observed from three different conc. of Zn and Cu solutions

Metal	Concentration (N)	Metal	Concentration (N)	EMF observed (E <sub>cell</sub> / V)
Zn/Zn <sup>2+</sup>	0.01 N	Cu/Cu <sup>2+</sup>	0.01 N	4-1.000
	0.02 N		0.02 N	7-1.085
	0.05 N		0.05 N	+1.086
			Average	+1.057

Results:

- (a). Standard electrode potential of Copper (E°) =  $\frac{+0.323}{}$  vs. SCE
- (c). EMF of the constructed Daniel cell = \_\_\_\_

Experimental Value	Actual Value	Percentage of error	Marks awarded
			11/2/12
			0/02/10
	-	1 Actual value	Actual Value