

Calculations:

A). Standardization of 0.1 N BaCl₂:

(N x V) of BaCl₂ solution = (N x V) of sodium sulphate

N of BaCl₂ solution = 0.02 N x 20 mL

Volume measured from Plot-1 (V₁) ; V₁ = 4.1 mL
= 0.097 N of BaCl₂ solution

B). Estimation of unknown sulphate:

(N x V) of irrigation water sample = (N x V) of BaCl₂ solution

N of irrigation water sample = N of BaCl₂ x Volume measured from Plot-2 (V₂) , V₂ = 4.5 mL
20 mL
= 0.022 N of irrigation water sample

Amount of sulphate present in 1L = Normality of irrigation water sample x Eq. wt. of SO₄²⁻ (48.03)

Amount of sulphate present in given sample solution = Strength of irrigation water sample x 48.03 x 100
1000
= 0.104 grams in 100 mL

Result: Amount of sulphate in given irrigation water sample = 0.104 grams.

Evaluation of Result:

Sample number	Experimental value	Actual Value	Percentage of error	Marks awarded